



Dyna-Form™ Crib Mattress

The **Dyna-Form Crib Mattress** has been specifically designed for use with newborn children. Its high frequency welded cover provides an excellent barrier from fluid ingress. This product can be made to measure in order to fit most crib & cot bed designs. Whilst using a special low odour foam, the application for use with infants makes this product particularly attractive to nursing staff. **Dyna-Form's** design quality comprises of a specialist high density CMHR foam surface, complete with a high frequency welded, multi-stretch, vapour permeable cover.

Features

Features	Benefits
① CMHR pressure reducing foam insert	Provides excellent comfort & pressure re-distribution properties
② Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear & friction forces, whilst assisting infection control policies
③ High frequency welded cover	Prevents fluid ingress & aids decontamination processes (can be laundered up to 80°C)
④ Non-slip base	Prevents slippage, maintaining patient safety
⑤ Fully replaceable parts	Ensures great value for money

Technical Specifications

Product Code	MAT/CRIB/SPECIAL
Risk Category	High Risk
Weight Limit	7 Stone / 44kg
Warranty	2 Years
Dimensions	Variable Sizes Available Upon Request
Foam Density	Beige – 36kg/m ³ Nominal Hardness – 125 Newtons
Fatigue Class	Very Severe
Fire Retardancy	BS7177 (Crib 5 Composite)
Mattress Weight	8kg
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)