



Dyna-Pad™ Leg Troughs

The **Dyna-Pad Leg Trough** range provides total positional care, with or without full ankle support, subject to product selection as below. The unique **Heel Free** trough immobilises the full leg from thigh to ankle ensuring maximum comfort, with a removed heel section for absolute pressure relief in an increasingly vulnerable area. Complete with a high frequency welded cover (and inspection zip for audit purposes) reduced shear and friction forces are achieved, whilst maintaining strict infection control standards.

Dyna-Pad™ Wedge

The **Dyna-Pad Wedge** provides knee support without the aid of an electric profiling bed. The foam wedge prevents the patient from slipping down the bed, whilst reducing shear & friction forces on the ankles and vulnerable heel area, therefore providing improved pressure re-distribution. Complete with an inspection zip for audit purposes, the multi-stretch, vapour permeable PU cover also reduces shear and friction forces.

Features

- ① Specialist CMHR foam cut with "heel free" design
- ② Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)
- Inspection zip
- ③ High frequency welded cover as standard

Benefits

- Maintains limb stability with complete pressure care, leaving the heel area free from pressure, shear & friction forces
- Reduces shear & friction forces, whilst assisting infection control policies
- Enhances audit procedures for infection control purposes
- Prevents fluid ingress & aids decontamination processes (can be laundered up to 95°C)

Technical Specifications

Product Code	1. LEG/HEEL/SIZE 2. LEG/SIZE 3. WEDGE/SIZE
Risk Category	High Risk
Warranty	1 Year
Dimensions	Variable Sizes Available Upon Request
Fire Retardancy	BSEN 597-1:1995 BSEN 597-2:1995
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)