



## Dyna-Pad™ Leg Troughs

The **Dyna-Pad Leg Trough** range provides total positional care, with or without full ankle support, subject to product selection as below. The unique **Heel Free** trough immobilises the full leg from thigh to ankle ensuring maximum comfort, with a removed heel section for absolute pressure relief in an increasingly vulnerable area. Complete with a high frequency welded cover (and inspection zip for audit purposes) reduced shear and friction forces are achieved, whilst maintaining strict infection control standards.

## Dyna-Pad™ Wedge

The **Dyna-Pad Wedge** provides knee support without the aid of an electric profiling bed. The foam wedge prevents the patient from slipping down the bed, whilst reducing shear & friction forces on the ankles and vulnerable heel area, therefore providing improved pressure re-distribution. Complete with an inspection zip for audit purposes, the multi-stretch, vapour permeable PU cover also reduces shear and friction forces.

### Features

- ① Specialist CMHR foam cut with "heel free" design
  - ② Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)
  - ③ High frequency welded cover as standard
- Inspection zip

### Benefits

- Maintains limb stability with complete pressure care, leaving the heel area free from pressure, shear & friction forces
- Reduces shear & friction forces, whilst assisting infection control policies
- Prevents fluid ingress & aids decontamination processes (can be laundered up to 80°C)
- Enhances audit procedures for infection control purposes

### Technical Specifications

<b>Product Code</b>	<b>1. LEG/HEEL/SIZE</b> <b>2. LEG/SIZE</b> <b>3. WEDGE/SIZE</b>
<b>Risk Category</b>	High Risk
<b>Warranty</b>	1 Year
<b>Dimensions</b>	Variable Sizes Available Upon Request
<b>Fire Retardancy</b>	BSEN 597-1:1995 BSEN 597-2:1995
<b>Hypochlorite Solution compatible</b>	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)