




Dyna-Tek™ Gel Cushion

The **Dyna-Tek Gel Cushion** has been designed using the latest polyurethane gel. By incorporating the most advanced raw materials during the manufacturing process, enhanced postural management and greater pressure re-distribution is achieved. The cooling gel surface assists any modern pressure ulcer management process by providing a minor reduction of the skin surface temperature. This product is specifically made for users considered to be at **“Very High Risk”** of pressure ulcer development. Equally at home in a ward, homecare or wheelchair environment, **Dyna-Tek Gel** gives advanced comfort and greater care to the user.

Features

- ① Centralised pommel design
- ② Complete gel surface
- ③ Finest quality CMHR foam
- ④ Multi-stretch, vapour permeable cover (impregnated with an anti fungal agent)
- ⑤ Non-slip base
- ⑥ Fully replaceable parts

Benefits

- Enhances postural management
- Provides excellent pressure re-distribution & important patient comfort
- Prevents patient bottoming out and aids patient transfer
- Reduces shear and friction forces, whilst assisting infection control policies (can be laundered up to 95°C)
- Aids moving & handling procedures, providing ultimate patient safety
- Ensures great value for money

Technical Specifications

Product Code	CUSH/GEL/46/46/7.5
Risk Category	Very High Risk
Weight Limit	25 Stone / 159kg
Warranty	2 Years
Dimensions	(Length) 46cm x (Width) 46cm x (Height) 7.5cm (18" x 18" x 3")
Foam Density	Blue Foam – 34kg/m ³ Nominal Hardness – 100 Newtons
Gel Type	Solid Polyurethane
Fatigue Class	Severe
Fire Retardancy	BS7176; 2007
Cushion Weight	3kg
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)