



## Dyna-Form<sup>®</sup> Wedge

The **Dyna-Form Wedge** provides knee support without the aid of an electric profiling bed. The foam wedge prevents the patient from slipping down the bed, whilst reducing shear & friction forces on the ankles and vulnerable heel area, therefore providing improved pressure re-distribution. Complete with an inspection zip for audit purposes, the multi-stretch, vapour permeable PU cover also reduces shear and friction forces.

### Features

①

High frequency welded cover as standard

Inspection zip

Specialist CME foam

Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)

### Benefits

Prevents fluid ingress & aids decontamination processes (can be laundered up to 95 °C)

Enhances audit procedures for infection control purposes

Maintains limb stability, free from pressure, shear & friction forces

Reduces shear & friction forces, whilst assisting infection control policies

### Technical Specifications

<b>Product Code</b>	<b>WEDGE/SIZE</b>
<b>Risk Category</b>	High Risk
<b>Warranty</b>	1 Year
<b>Dimensions</b>	Variable Sizes Available Upon Request
<b>Fire Retardancy</b>	BSEN 597-1:1995 BSEN 597-2:1995
<b>Colour</b>	Blue
<b>Hypochlorite Solution compatible</b>	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)