




Dyna-Flex™ Bari Cushion

The **Dyna-Flex™ Bariatric Cushion** incorporates the most exclusive combination of visco elastic & CME foam technology, providing not only outstanding pressure re-distribution, but great patient comfort. The temperature sensitive memory foam softens at patient interface, thus giving even weight distribution, whilst the CME underside prevents the foam from bottoming out. Suitable for larger patients considered to be at **“High Risk”** of pressure ulcer development, the **Dyna-Flex** gives nursing security without the patient challenging the issue of comfort. Complete with a well designed cover and non-slip base, this product is ideal for use in conjunction with any modern wheelchair or general seating environment.

Features

- | | | |
|---|--|--|
| ① | Temperature sensitive memory foam | Moulds to patient shape, providing excellent comfort & pressure re-distribution properties |
| ② | Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent) | Reduces shear & friction forces, whilst assisting infection control policies |
| ③ | High frequency welded cover | Prevents fluid ingress and aids decontamination processes (can be laundered up to 80°C) |
| ④ | Non-slip base | Aids moving & handling procedures, providing ultimate patient safety |
| ⑤ | Fully replaceable parts | Ensures great value for money |
| ⑥ | Available with sag compensator | Suitable for wheelchair users |

Benefits

Technical Specifications

Product Code	CUSH/FLEX/BARI/SIZE
Risk Category	High Risk
Weight Limit	50 Stone / 318kg
Warranty	1 Year
Dimensions	Made to Measure / Variable Sizes Available Upon Request
Foam Density	Visco 55 - 60kg/m3 – Nominal Hardness Soft to Medium CME Base 38 - 40kg/m3 – Nominal Hardness 200 Newtons
Fatigue Class	Severe
Fire Retardancy	BS7176
Weight	1.2kg
Colour	Wine
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)