



Dyna-Tek™ Pad Cushion

The **Dyna-Tek Pad Cushion** has been created for users considered to be **“At Risk”** of pressure ulcer development. Dyna-Tek’ design quality comprises of a specialist “castellated cut” CME foam surface, ensuring the reduction of heat & moisture build up whilst achieving maximum pressure re-distribution. Complete with a sewn, multi-stretch, vapour permeable cover and non-slip base, this economical product is ideal for use in a Wheelchair, Hospital or Homecare environment.

Features

①	CME castellated foam cushion (with tear duct cut)	Provides excellent comfort and pressure re-distribution properties
②	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear & friction forces, whilst assisting infection control policies
③	Sewn or welded cover (NB: welded cover option carries a small price premium)	Prevents fluid ingress and aids decontamination processes (can be laundered up to 95°C)
④	Non-slip base (NB: carries a price premium)	Prevents slippage, encouraging patient safety
⑤	Available with sag compensator	Suitable for wheelchair users

Benefits

Technical Specifications

Product Code	CUSH/PAD/43/43/10
Risk Category	At Risk
Weight Limit	20 Stone / 127kg
Warranty	2 Years Cover, 4 Years Foam
Dimensions	(Length) 43cm x (Width) 43cm x (Height) 10cm (17" x 17" x 4") Variable Sizes Available Upon Request
Foam Density	34–36kg/m ³ Nominal Hardness – 125 Newtons
Fatigue Class	Very Severe
Fire Retardancy	BS7176
Weight	1kg
Colour	Light Blue
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)