




# Dyna-Tek™ Profile

The **Dyna-Tek Profile** Cushion is our entry level product, providing surprisingly good postural management, along with great pressure re-distribution. This product is specifically made for users considered to be **“At Risk”** of pressure ulcer development. Designed mainly for use in a wheelchair but equally applicable to that of a ward or homecare environment, **Dyna-Tek Profile** provides the ultimate combination of pressure care and patient comfort at a cost effective price.

## Features

- |   |   |  |
|---|---|--|
| ① | Unique “ischial tuberosity well” design                                       | Creates exceptional postural management  |
| ② | Sculpted thigh support foam surface   | Provides excellent pressure re-distribution & important patient comfort                                    |
| ③ | Finest quality CME foam   | Prevents patient bottoming out and aids patient transfer   |
| ④ | Multi-stretch, vapour permeable cover (impregnated with an anti fungal agent) | Reduces shear & friction forces, whilst assisting infection control policies (can be laundered up to 95°C) |
| ⑤ | Non-slip base   | Aids moving & handling procedures, providing ultimate patient safety                                       |
| ⑥ | Fully replaceable parts   | Ensures great value for money  |

## Benefits

## Technical Specifications

<b>Product Code</b>	<b>CUSH/PRO/43/43/5</b>
<b>Risk Category</b>	At Risk
<b>Weight Limit</b>	18 Stone / 114kg
<b>Warranty</b>	1 Year
<b>Dimensions</b>	(Length) 43cm x (Width) 43cm x (Height) 5cm (17” x 17” x 2”) Variable Sizes Available Upon Request
<b>Foam Density</b>	CME Base – 38–40kg/m <sup>3</sup> Nominal Hardness – 200 Newtons
<b>Fatigue Class</b>	Severe
<b>Fire Retardancy</b>	BS7176; 2007
<b>Weight</b>	1kg
<b>Colour</b>	Light Blue
<b>Hypochlorite Solution compatible</b>	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)