



The **Dyna-Flex Cushion** incorporates the most exclusive combination of visco elastic & CME foam technology, providing not only outstanding pressure re-distribution, but great patient comfort. The temperature sensitive memory foam softens at patient interface, thus giving even weight distribution, whilst the CME underside prevents the foam from bottoming out. Suitable for patients considered to be at **“High Risk”** of pressure ulcer development, **Dyna-Flex** gives nursing security without the patient challenging the issue of comfort. Complete with a well designed cover and non-slip base, this product is ideal for use in all hospital, nursing or homecare environments.

Features

①	Temperature sensitive memory foam	Moulds to patient shape, providing excellent comfort & pressure re-distribution properties
②	Multi-stretch, vapour permeable fabric (impregnated with an anti fungal agent)	Reduces shear & friction forces, whilst assisting infection control policies
③	High frequency welded cover	Prevents fluid ingress and aids decontamination processes (can be laundered up to 95°C)
④	Non-slip base	Prevents slippage, encouraging patient safety
⑤	Fully replaceable parts	Ensures great value for money
⑥	Available with sag compensator	Suitable for wheelchair users

Benefits

Technical Specifications

Product Code	CUSH/FLEX/43/43/7.5
Risk Category	High Risk
Weight Limit	35 Stone / 222kg
Warranty	2 Years
Dimensions	(Length) 43cm x (Width) 43cm x (Height) 7.5cm (17" x 17" x 3") Variable Sizes Available Upon Request
Foam Density	Visco – 55–60kg/m ³ Nominal Hardness – Soft to Medium CME Base – 38–40kg/m ³ Nominal Hardness – 200 Newtons
Fatigue Class	Severe
Fire Retardancy	BS7176; 2007
Weight	2kg
Colour	Wine
Hypochlorite Solution compatible	1,000 ppm diluted or 0.1% (or to be used in conjunction with local infection control guidelines)